

## Non-commercial Plasmid Mini Preparations

Ilsa León, Llinás Lab 2007

### Reagents

- **Alkaline Lysis Buffer:** 0.1 M Tris-HCl pH 8.0; 0.1 M NaOH; 0.5% SDS; 1 mM EDTA; add RNase A to 50 µg/ml prior to use from a 10 mg/ml stock stored in -20°C  
**For 10 ml:** 1 ml 1 M Tris-HCl, pH 8.0; 1 ml 1 M NaOH; 20 µl 0.5 M EDTA; 500 µl 10% SDS. Add 5 µl of a 10 mg/ml RNase A stock for every 1 ml of Alkaline Lysis Buffer needed in the below protocol.
- 3 M Na-Acetate pH 5.3
- Isopropanol (-20°C)
- Tris-EDTA (TE) buffer
- Phenol:Chloroform (1:1 v/v, pH of phenol must be between 7.8-8.0)
- 100% ethanol
- 70% ethanol
- LB

### A. Preparation of Cells

1. Start a 25 ml O/N culture with the proper antibiotic.
2. Take 1.5 ml of culture and pellet the cells by centrifuging for 60 seconds in a microcentrifuge. Remove remaining media by aspiration, leaving the bacterial pellet as dry as possible. Repeat 2 more times.

### B. Lysis of Cells

3. Add 400 µl alkaline lysis buffer and loosen the pellet with vigorous vortexing to lyse the cells. (Run your Eppendorf tubes along a microfuge rack several times to assist you. It is very important that you resuspend your pellet completely). Place on ice for 5 min.
4. Add 200 µl 3 M Na-Acetate pH 5.3 and mix well by briefly vortexing.
5. Centrifuge sample at maximum speed in a tabletop centrifuge for 5-10 min.

### C. Recovery of Plasmid DNA

6. Transfer supernatant to a fresh tube and precipitate by adding 550 µl cold isopropanol. Mix well by vortexing. Incubate for 10 min on ice and centrifuge for 30 min at 4°C, max speed (or incubate in -20 °C indefinitely).
7. Dry the pellet briefly and resuspend in 500 µl TE or water.

### D. Phenol:Chloroform Extraction

8. Add an equal volume of phenol:chloroform.
9. Mix by vortexing until an emulsion forms.
10. Centrifuge the emulsion at maximum speed for 2 minutes.
11. Transfer the aqueous upper layer to a fresh tube without pipetting the interface or organic phase.
12. Repeat steps 8-11 with the extracted aqueous layer.

### Materials

- 50 ml conical tubes
- 1.5 ml tubes
- Allegra® 25R and TA10.250 rotor plus yellow adapters for 15 ml conical tubes
- Microcentrifuge in 4°C
- Microcentrifuge at room temperature
- Vortex
- Disposable pipette tip attached to vacuum line and side armed flask
- Laminar air flow cabinet

### **E. Ethanol Precipitation**

13. Add 3 M Na-Acetate pH 5.3 at 1/10 the volume of the aqueous phase.
14. Mix the solution well.
15. Add 3 volumes of (ice-cold) 100% ethanol and mix the solution well.
16. Incubate the ethanol solution on ice for 15-30 minutes to allow the DNA to precipitate (this solution may be left indefinitely at 0°C or at -20°C).
17. Pellet DNA by centrifugation at max speed in 4°C for 30 min. Carefully remove the supernatant by aspiration.
18. Wash the pellet with 500 µl of 70% ethanol. Centrifuge the sample at max speed for 2 minutes at 4°C. Carefully remove the supernatant as above.
19. Perform the following steps in a laminar air flow cabinet. Dry the pellet at room temperature until the last traces of fluid have been evaporated. The pellet represents approximately 40-60 µg DNA.
20. Resuspend the pellet in appropriate amount of cytomix and store in 4°C until ready to use.
21. Spec 1:100 dilution checking OD<sub>260</sub> and OD<sub>260/280</sub>.