

# Polymerase Chain Reaction for *Plasmodium falciparum*

Ilsa León, Llinás Lab 2007

## Reagents

- Primers 100 picomol/ $\mu\text{l}$ <sup>1</sup>
- 25 mM dNTPs (1:2:2:1, G:A:T:C)<sup>2</sup>
- 25 mM MgCl<sub>2</sub>
- 10X Taq Buffer (500 mM KCl, 100 mM Tris-HCl, pH 8.3)
- 0.1-1 nanogram gDNA
- Taq or other DNA polymerase
- Milli Q water (or ddH<sub>2</sub>O)
- Agarose
- Tris-Acetate EDTA (TAE) or other electrophoresis buffer
- DNA ladder
- Ethidium Bromide

## Materials

- 0.2 ml PCR tubes
- Thermal cycler
- Microwave
- Electrophoresis apparatus
- Gel imaging apparatus

1. When designing primers stay within the lower end of the range for T<sub>m</sub> [60-75°C] keeping both the forward and reverse primers within 5°C of one another. Measure this temperature only for the base pairs that actually will be annealing during the first annealing step of the PCR program.
2. dNTPs are usually sold separately as 100 mM. Use the above proportions and bring the final total concentration of dNTPs to 25 mM using Milli Q. One example is the following mix: 20  $\mu\text{l}$  A, 20  $\mu\text{l}$  T, 10  $\mu\text{l}$  G, 10  $\mu\text{l}$  C, 180  $\mu\text{l}$  Milli Q.

## PCR Setup

0.5  $\mu\text{l}$  Forward Primer  
0.5  $\mu\text{l}$  Reverse Primer  
0.5  $\mu\text{l}$  dNTPs  
2.5  $\mu\text{l}$  MgCl<sub>2</sub>  
5.0  $\mu\text{l}$  Taq Buffer  
5.0  $\mu\text{l}$  0.1 ng/ $\mu\text{l}$  gDNA  
35.0  $\mu\text{l}$  Milli Q  
1.0  $\mu\text{l}$  Taq

## PCR Program A

1. 92°C 1 min  
2. 92°C 30 sec  
3. 54°C 30 sec  
4. 62°C 1 min 30 sec  
5. Go to step 2, 29 times  
6. 62°C 4 min  
7. 4°C FOREVER

## PCR Program B

1. 94°C 1 min  
2. 94°C 30 sec  
3. 65°C 30 sec  
4. 72°C 30 sec  
5. Go to step 2, 29 times  
6. 72°C 4 min  
7. 4°C FOREVER

\*Run 4  $\mu\text{l}$  of PCR reactions on an agarose gel using a DNA Standard to compare fragment sizes, ethidium bromide to stain DNA, and an UV gel imager to visualize the bands.