

**\*\* STERILE FILTER ALL BUFFERS TO AVOID DUST ON ARRAYS\*\***

	Final Concentration	10x PBS	Milli-Q H <sub>2</sub> O	Other
<b>1x PBS</b>	---	100ml	900ml	---
<b>Solution A</b>	0.01% TX-100 / PBS	100ml	899ml	1ml 10% TX-100 (in H <sub>2</sub> O)
<b>Solution B</b>	0.1% Tween-20 / PBS	100ml	895ml	5ml 20% Tween-20 (in H <sub>2</sub> O)
<b>Solution C</b>	0.5% Tween-20 / PBS	100ml	875ml	25ml 20% Tween-20 (in H <sub>2</sub> O)
<b>Solution D</b>	0.05% Tween-20 / PBS	100ml	897.5ml	2.5ml 20% Tween-20 (in H <sub>2</sub> O)

1. Dissolve milk solutions on shaker for 2 hours (200rpm), syringe filter 0.45um and keep on ice
  - a. 2% milk: 0.1g powdered milk in 5ml 1X PBS
  - b. 4% milk: 0.1g powdered milk in 2.5ml 1X PBS
2. Prepare Protein Binding Mixture:
  - a. Thaw protein and reagents on ice, label tube with protein and array number

**PROTEIN-BINDING MIXTURE:**

18.6ul Milli-Q H<sub>2</sub>O (sterile)  
 15.7ul 4X PBS (sterile)  
 62.5ul 4% milk in 1X PBS (sterile)  
 0.7ul Salmon testes DNA (~53 ng/ul stock)  
 2.5ul 10 mg/ml BSA (0.2 ug/ul final conc.)  
 25ul Protein (in elution buffer with 25% glycerol)  
 125ul

- b. Flick tube to mix and spin briefly to get rid of bubbles
  - c. Incubate mixture on ice for 1 hour
    - i. *Note - leftover protein cannot be re-frozen or used again*
    - ii. *If working with multiple samples, make a master-mix of all reagents except protein*
3. Block for non-specific binding:
  - a. Pre-wet double-stranded array in **Solution A** in Coplin jar (max. 2 arrays per jar)
    - i. Shake on platform at 125rpm for 5 min (cover with foil)
  - b. Add 800ul 2% milk to single-chamber gasket slide positioned in steel hybe chamber
    - i. Use pipette tip to spread the milk evenly over the well without touching gasket
  - c. Fill a staining dish with 1X PBS (use throughout the entire protocol, keep covered)
  - d. Rinse array briefly in PBS and pull out of solution slowly to remove droplets
  - e. Lower array ("Agilent" side facing down) onto the gasket and close hybe chamber
    - i. Be careful to avoid bubbles and avoid wetting the top of the gasket
  - f. Block array in 2% milk for 1 hour in the dark (cover the chamber with foil or ice bucket)
4. Wash off blocking agent:
  - a. Fill a staining dish and Coplin jar with **Solution B**
  - b. Remove gasket/slide sandwich, taking care not to break the seal
  - c. Disassemble while submerged in Solution B using plastic tweezers
  - d. Shake slide underwater and quickly transfer to Coplin jar (shake @ 125rpm for 5 min)
    - i. *Rinse staining dish and gasket slide with water, spin slide dry and return to box*
  - e. Repeat Coplin jar wash step using **Solution A** (shake 2 min)

5. Bind Protein:

- a. Prepare steel hybe chamber with 4-well gasket slide, note orientation:
- b. Pipette 125ul protein binding mixture into each well
  - i. *Note - it is best to use 2 non-adjacent wells per slide (i.e. 1 + 3 OR 2 + 4) to minimizes the damage done if one of the gaskets leaks*
- c. Remove slide from Coplin jar, rinse briefly in PBS as before and remove slowly
- d. Lower array ("Agilent" side facing down) onto the gasket slide and close hybe chamber
- e. Incubate protein binding mixture on array for EXACTLY 1 hour in the dark



6. During the protein binding step, prepare fluorescent antibody mixture (will need 125ul/array):

**FLUORESCENT ANTIBODY MIXTURE:**

# of Arrays	2% Milk	Alexa488-conjugated anti-GST antibody (4°C)	Total Vol.
1	125.8ul	3.2ul	129ul
2	248.6ul	6.4ul	255ul
3	370.5ul	9.5ul	380ul
4	492.4ul	12.6ul	505ul

- a. Wrap tube in foil and incubate antibody mixture on ice for 1 hour
  - b. *This is a 1:40 dilution of primary antibody (Molecular Probes #A-11131, 2 mg/ml stock).*
7. Wash off protein binding mixture:

- a. Fill a staining dish and Coplin jar with **Solution C**
- b. Remove gasket/slide sandwich, taking care not to break the seal
- c. Disassemble while submerged in Solution C using plastic tweezers
- d. Shake slide underwater and quickly transfer to Coplin jar (shake @ 125rpm for 3 min)
- e. Repeat Coplin jar wash step using **Solution A** (shake 2 min)
  - i. *Wash 4-well gasket slide with water and spin slide dry*

8. Bind Fluorescent Antibody:

- a. Prepare hybe chamber with dried 4-well gasket slide and 125ul antibody mixture
- b. Remove slide from Coplin jar, rinse briefly in PBS as before, removing slide from PBS slowly
- c. Lower array ("Agilent" side facing down) onto the gasket and close hybe chamber
- a. Incubate fluorescent antibody mixture on array for 1 hour in the dark

9. Wash off fluorescent antibody mixture:

- a. Fill a staining dish and two Coplin jars with **Solution D**
- b. Remove gasket/slide sandwich, taking care not to break the seal
- c. Disassemble while submerged in Solution D using plastic tweezers
- d. Shake slide underwater and quickly transfer to Coplin jar (shake @ 125rpm for 3 min)
- e. Repeat Coplin jar wash step using fresh **Solution D** (shake 3 min)
- f. Repeat Coplin jar wash step using **1x PBS** (shake 2 min)
- g. Rinse briefly in PBS staining dish as before, pulling up slowly to remove liquid droplets
  - i. Spin array dry if droplets remain

10. Store array slide in a slide box in the dark until ready to scan

*Note: PBM replicates should be run on different versions – same probes, different arrangement*

## Scanning PBM Arrays

1. Scan arrays using an Axon 4300A scanner or similar (requires 488nm blue laser)
  - a. Wavelength #1: 488nm
  - b. Pixel Size: 5um
  - c. PMT Gain: 350
  - d. Power (%): 100
  - e. Filter: Standard Blue
  - f. Lines to Average: 1
2. Scan full slide at low resolution to preview using “>” button
3. Select single array and scan at high resolution using “>” button
  - a. Use “!” button to auto-contrast image
  - b. Zoom in to check Simple Line Profile
    - i. Click on Measuring Tools (pencil icon) and select Simple Line Profile
    - ii. Draw a line across a row of spots with at least one bright (saturated) spot
    - iii. Graph should show a large peak representing bright spot
  - c. Once you get saturated white spots (around 60,000), Save Image – Selected Scan Area
    - i. [20130325\_v11#417\_ProteinName\_Alexa488\_PMT350.tif]
  - d. Repeat scan at lower PMT Gain using the Simple Line Profile readout as a guide
    - i. Highest scan peak = 60,000
    - ii. Middle scan peak = 40,000
    - iii. Lowest scan peak = 5,000-10,000
4. Gene Pix Pro software:
  - a. Open lowest PMT scan and import settings file from double-stranded scan (.gps file)
  - b. Align the grid and flag spots for removal
  - c. Save Settings As [20130325\_v11#417\_ProteinName\_Alexa488.gps]
  - d. Click “Analyze” to generate a results file for each of the 4 arrays
  - e. Repeat above steps for the middle and high-PMT scans

**For more information:** Berger MF and Bulyk ML. Universal protein-binding microarrays for the comprehensive characterization of the DNA-binding specificities of transcription factors. *Nature Protocols*. (2009) 4(3):393-411

### Folder: PBM FILES > v11#417

20130125\_v11#417\_Cy3.gps

*Double stranding settings file (grid alignment and flagged spots)*

20130125\_v11#417\_Cy3\_PMT350.tif

*Double stranding high-res scan at 350 PMT*

20130125\_v11#417\_Cy3\_PMT350\_1-4.gpr

20130125\_v11#417\_Cy3\_PMT350\_2-4.gpr

20130125\_v11#417\_Cy3\_PMT350\_3-4.gpr

*Double-stranding fluorescent values for each spot including info on what was flagged*

20130125\_v11#417\_Cy3\_PMT350\_4-4.gpr

### Folder: PBM FILES > v11#417 > v11#417-1 ProteinA

20130125\_v11#417\_ProteinA\_Alexa488.gps

*Alexa488 settings file (grid alignment and flagged spots)*

20130325\_v11#417\_ProteinA\_Alexa488\_PMT350.tif

*Alexa488 high-res scan at 350 PMT*

20130325\_v11#417\_ProteinA\_Alexa488\_PMT300.tif

*Alexa488 high-res scan at 300 PMT*

20130325\_v11#417\_ProteinA\_Alexa488\_PMT250.tif

*Alexa488 high-res scan at 250 PMT*

20130325\_v11#417\_ProteinA\_Alexa488\_PMT350\_1-4.gpr

*Alexa488 quantitative data for array 1*

20130325\_v11#417\_ProteinA\_Alexa488\_PMT350\_2-4.gpr

*Alexa488 quantitative data for array 2 (junk)*

20130325\_v11#417\_ProteinA\_Alexa488\_PMT350\_3-4.gpr

*Alexa488 quantitative data for array 3 (junk)*

20130325\_v11#417\_ProteinA\_Alexa488\_PMT350\_4-4.gpr

*Alexa488 quantitative data for array 4 (junk)*

20130325\_v11#417\_ProteinA\_Alexa488\_PMT300\_1-4.gpr

*Alexa488 quantitative data for array 1*

20130325\_v11#417\_ProteinA\_Alexa488\_PMT300\_2-4.gpr

*Alexa488 quantitative data for array 2 (junk)*

20130325\_v11#417\_ProteinA\_Alexa488\_PMT300\_3-4.gpr

*Alexa488 quantitative data for array 3 (junk)*

20130325\_v11#417\_ProteinA\_Alexa488\_PMT300\_4-4.gpr

*Alexa488 quantitative data for array 4 (junk)*

20130325\_v11#417\_ProteinA\_Alexa488\_PMT250\_1-4.gpr

*Alexa488 quantitative data for array 1*

20130325\_v11#417\_ProteinA\_Alexa488\_PMT250\_2-4.gpr

*Alexa488 quantitative data for array 2 (junk)*

20130325\_v11#417\_ProteinA\_Alexa488\_PMT250\_3-4.gpr

*Alexa488 quantitative data for array 3 (junk)*

20130325\_v11#417\_ProteinA\_Alexa488\_PMT250\_4-4.gpr

*Alexa488 quantitative data for array 4 (junk)*